

Wheelchair Lift Vehicle Lift



This easy to use product lifts and stores your unoccupied manual folding wheelchair or travel scooter* inside your vehicle at the touch of a button. Simply connect the docking device, push a button, and let your Wheelchair Lift do the heavy lifting for you!

Model AWL-150 Highlights:

- Power Raise/Lower and Manual Rotate Operation - Makes lifting and storing your manual folding wheelchair or travel scooter inside your vehicle virtually effortless (150 lb/68 kg lifting capacity)
- Lift folds down when not in use
- Easy, one button operation
- Docking device for your mobility device is included at no charge
- Safe/Reliable - lifecycle and strength tested beyond industry standards
- No structural modifications to your vehicle - can be easily removed for vehicle resale
- Can reinstall in your next vehicle
- Works with minivans, SUVs, CUVs, full-size vans and sedans
- Hundreds of custom, vehicle-specific mounting kits are available for an extra clean and precise installation
- Bruno Turning Automotive Seating™ (TAS™) complements your Wheelchair Lift for a total mobility solution

Important Note: Confirm all applications

Rochester, NY
585.385.2556
951 Panorama Trail
Rochester, NY 14625

Syracuse, NY
315.510.4272
1241 Erie Blvd. West
Syracuse, NY 13204